



**Improving the  
lives of people  
suffering from  
Dry Eye  
Syndrome**

**Treat The  
Cause of Dry  
Eye**



[TearRestore.com](http://TearRestore.com)



**TearRestore**

**First Open-Eye  
Warm Compress**

**Unique Convenience, Effective  
Heat, Natural Meibomian Gland  
Expression**

# Dry Eye Syndrome

Tears are necessary for maintaining the ocular surface and providing clear vision. Dry Eye Disease is a condition in which a person is unable to produce an adequate quality or quantity of tears to lubricate and nourish the eye. The symptoms include burning, itching, eye fatigued, redness, light sensitivity, and blurred vision.

# Meibomian Gland Dysfunction (MGD)

MGD is the primary cause of Dry Eye Syndrome and occurs when the oil glands within the eyelids, Meibomian glands, are not working properly and become clogged. The treatment for MGD is a daily warm compress. This treatment applies heat to the eyelids for 10 minutes, and results in the softening of the clogged oils in the eyelids. This allows the oils to flow more freely and helps prevent the tears from evaporating from the surface of the eye.

# TearRestore Solution

TearRestore's open-eye warm compress delivers an effective, user friendly at-home treatment for MGD. Developed with assistance from industry leading physicians, TearRestore delivers therapeutic heat where it is needed while permitting the user to see throughout treatment.



# How it Works

The TearRestore mask has specially designed reusable heat packs that target the Meibomian glands and supply therapeutic heat for a minimum of ten minutes every use. Additionally, TearRestore is the only open-eye warm compress which allows you to be active during your daily warm compress! With TearRestore, it is easy to find time to treat your eyes everyday.



## Step I "Click" the Metal Disc

Within the heat pack there is a small metal disc. Press the center of the metal disc with your thumb or forefinger; this will activate the pack and you will begin to feel the heat.

## Step II Massage the Heat Pack and Attach the Heat Pack to the Mask

Massage the heat pack after it is fully activated to evenly spread the heat. Next, match the eye holes of the heat pack with the mask and hook the edges of the heat pack eye holes onto the mask.

## Step III Adjust For Comfort

Place the TearRestore mask and activated heat pack over your head and adjust to be able to see through the eye holes. Tighten the adjustable straps until comfortable then relax and/or continue your day!

# Reactivation

The TearRestore heat packs are easily reactivated in boiling water or the TearRestore "one-click kettle". Simply boil enough water to submerge the heat packs and place the heat packs in the boiling water for approximately 5 minutes. Once the heat packs have returned to a fully liquid state, remove them from the boiling water and let cool.



# Improve Patient Outcomes

The TearRestore mask offers patients the convenient treatment they want while significantly improving the signs and symptoms of dry eye syndrome.

# Consistent Heat

TearRestore supplies therapeutic heat (>40C) for a minimum of ten minutes every use. This ensures patients and doctors an effective treatment every time.

# Targeted Treatment


TearRestore's design utilizes anatomical norms to target the oil glands of the eyelid while avoiding the globe. This permits the user to blink throughout treatment, resulting in natural oil gland expression throughout use.

↑ 29% INCREASE IN OBJECTIVE LIPID LAYER THICKNESS


↑ 43% IMPROVEMENT IN SUBJECTIVE SPEED SCORES

\*BASED ON PILOT STUDY DATA

TearRestore Heat Pack Temperature



Thermal Image Post TearRestore



Temperature →

