



Suggested Vitreoretinal Set



6-676 DK Squeeze Handle for VR Instrument Heads

- Suitable for all Duckworth & Kent Vitreoretinal Instrument Heads
- Overall Length 92 mm (without head)
- Squeeze action handle



Scissors



Straight Scissors

23G

- Sharp tipped straight micro scissors with 2.2mm long blade
- Straight 23 gauge tube, length 31.3mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 144mm

1-841N



23 Gauge



Straight Blunt Scissors

23G

- Blunt tipped straight micro scissors with 1.9mm long blade
- Straight 23 gauge tube, length 31.3mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 144mm

1-841BN



23 Gauge



Curved Scissors

NEW tip design 23G

- Sharp tipped curved micro scissors with 2.2mm long blade
- Straight 23 gauge tube, length 31.3mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 144mm

1-842N



23 Gauge

Forceps



End Gripping Forceps
ILM Peeling

23G

- Fine end gripping tips, tip length 3mm
- Straight 23 gauge tube, length 31mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 142mm

2-877N



23 Gauge



Delicate Serrated
Forceps

23G

- Delicate serrated jaws with blunt tip, tip length 3.5mm
- Straight 23 gauge tube, length 31mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 142mm

2-878-1N



23 Gauge



Asymmetrical Forceps

23G

- 2mm jaws, 0.7mm platforms, 55° angled pointed picking tip
- Straight 23 gauge tube, length 31mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 142mm

2-879N



23 Gauge

Fixed Heads


23 Gauge Forceps & Scissors

The fixed head range of Vitreoretinal forceps and scissors are a one piece devices featuring the same tip and handle design as the interchangeable head.

All fixed head vitreoretinal forceps and scissors are supplied with a protective guard which can be sterilised.




Scissors



Straight Scissors

23G

- Sharp tipped straight micro scissors with 2.2mm long blade
- Straight 23 gauge tube, length 32.5mm
- Squeeze action handle closes scissors
- Round squeeze handle, length 138mm



1-841


23 Gauge



Straight Blunt Scissors

23G

- Blunt tipped straight micro scissors with 1.75mm long blade
- Straight 23 gauge tube, length 32.5mm
- Squeeze action handle closes scissors
- Round squeeze handle, length 138mm



1-841B


23 Gauge



Curved Scissors

23G

- Sharp tipped curved micro scissors with 2.5mm long blade
- Straight 23 gauge tube, length 32.5mm
- Squeeze action handle closes scissors
- Round squeeze handle, length 139mm



1-842

23 Gauge

Duckworth & Kent - Titanium

Forceps




End Gripping Forceps
ILM Peeling

23G

- Fine end gripping tips, tip length 3mm
- Straight 23 gauge tube, tube length 32.5mm
- Squeeze handle activates both jaws
- Round squeeze handle, overall length 140mm

2-877




23 Gauge

Delicate Serrated
Forceps

23G

- Delicate serrated jaws with blunt tip, tip length 3mm
- Straight 23 gauge tube, tube length 32.5mm
- Squeeze handle activates both jaws
- Round squeeze handle, overall length 139mm

2-878-1




23 Gauge

Asymmetrical Forceps

23G


- Asymmetrical tip design with fine platforms
- Straight 23 gauge tube, tube length 32.5mm
- Squeeze handle activates both jaws
- Round squeeze handle, overall length 138mm

2-879



23 Gauge


6-800 *Barrett LeClip Utility Clamp*



- 14mm serrated cross action jaws
- Length 82mm
- Distinctive identification labelling




6-805 *Bulldog Clip*



- 8.5mm serrated cross action jaws
- Length 46mm




9-650 *DK Castroviejo Style Marking Caliper*



- 0.25mm x 0.25mm delicate marking tips
- Marks 0mm to 15mm in 0.25 increments
- Marks from centre of tips
- Adjustable thumb screw
- Standard caliper handle




2-160 *Troutman Superior Rectus Forceps*



- 0.5mm, 1 x 2 teeth
- 45° angled shafts, tip to angle length 10mm
- Flat handle, length 113mm




2-118E *Bonn Suturing Forceps (0.3mm)*



- 0.3mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- Flat handle, length 90mm
- Design registration number 004383396-0001




2-132E Colibri Toothed Forceps (0.12mm)




- 0.12mm, 1 x 2 teeth, tip length 2mm
- 6mm tying platforms
- Colibri style shafts
- Flat handle, length 90mm
- Design registration number 004383396-0001




2-500E McPherson Tying Forceps



- 6mm tying platforms
- Straight shafts
- Flat handle, length 90mm
- Design registration number 004383396-0001




1-118 DK Castroviejo Miniature Corneal Scissors (curved)



- Blunt tips, curved blades
- Universal cutting blades
- Cut length 6mm, tip to pivot length 11mm
- Flat handle, length 115mm




1-510 Westcott Style Tenotomy Scissors (curved)



- Blunt tips, curved blades
- Cut length 10mm, blade length 20mm
- Flat handle, length 115mm



4-600 Pearce Retractable Diamond Knife, 45° Double Edge



W=1.0T=0.17L=3.5

- 1mm wide 45° double edge diamond blade
- Second position 10mm shaft extension
- Retractable long handle, Length 124mm



6-912 Giunchiglia Membrane Peeling Spatula, 23 Gauge



- Ultra fine textured tip
- 0.6mm wide, blunt tip curved 45°



- Straight 23 gauge shaft, length 32mm
- Round handle, overall length 128mm

6-912-1 Membrane Peeling Spatula, 23 Gauge




- Fine textured tip
- 0.6mm wide, blunt tip curved 45°





- Straight 23 gauge shaft, length 32mm
- Round handle, overall length 128mm

9-556 Kratz Barraquer Speculum

- 14.5mm open blades
- Angled to rest temporally
- Single piece construction




6-635-2 DK Schocket Double Ended Scleral Depressor



- 2.5mm round and 4.5mm cylinder tips
- Round handle, length 135mm




6-625 Jameson Muscle Hook



- 2mm bulbous tip, tip to curve length 8mm
- Round handle, length 136mm




3-302 DK Barraquer Needle Holder (curved with lock)



- 9mm delicate curved jaws
- Round handle with lock, hinge spring
- Length 109mm



3-203 DK Needle Holder (curved)






- 9mm medium curved jaws
- Round handle, tag spring
- Length 110mm




6-675 Protective Cleaning Cover

- Suitable for all Duckworth & Kent Vitreoretinal Instrument Heads
- Protects the Vitreoretinal detachable heads
- Standard luer fitting for flushing attachment



T7020 Sterilising Tray (suitable for 20 instruments)

- External dimensions: 395mm x 266mm x 25mm
- Suitable for 20 instruments
- One silicone mat



www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.co.uk



© April 2011 Duckworth & Kent
21.05.18

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeon's names or by the category of the instrument.

at the Leading Edge