


PalmScan™



Solutions That Fit

Ultra Mobile B-Scan

Introducing MMD's portable and most versatile ultrasound device. Turns any PC* into a B-Scan with high resolution imaging.





Features:

- **Easy To Use**
Custom settings for multiple applications
- **Patient Database**
Virtually unlimited
- **PC Adaptable-Travels With You***
Same probe, multiple PCs
- **Adjustable Pulse Power**
Improves penetration capability
- **High Resolution Imaging**
2,000 samples per line
- **Dual Calipers**
Easily measures structures within the eye
- **Printing Capability**
Simply print captured images to your existing office printers
- **Single Frame and Cine Modes**
Store movies or JPEG images with adjustable size cine buffers
- **Zoom Capability**
Directly zooms into the structure within the eye
- **Pulse Frequency 12MHz/15MHz**
3-10 cm penetration depth



* requires USB 2.0 interface

Micro Medical Devices, Inc.
(866) 730-0663 (OMMD)
www.micromedinc.com
info@micromedinc.com



Notes

Niente da segnalare

B-Scan Parameters

Caliper A		Caliper B		AMode Vector		AMode Gain
0.00		0.00		127		158
Depth	Frequency	Intensity	Contrast	Far Gain	Mid Gain	Near Gain
5 cm	15 MHz	44	80	22	87	10

Probe name: OP 12 MHz, ID: 06861