

blink[®] medical

by corzamedical

ophthalmic single use instruments and packs



www.blinkmedical.com

Welcome

Blink Medical is one of the world's largest manufacturers of quality, sterile, single use surgical instruments and bespoke instrument sets and packs. We supply instruments from our UK headquarters, to both NHS and private sector hospitals for use in either theatres or out patient departments. Through our exclusive distributor network, we export to over 20 countries worldwide including Australia, New Zealand and throughout continental Europe.

Offering the highest quality instruments, we provide products for Ophthalmic, Dermatology, Plastic Surgery, Cosmetic Surgery and ENT departments. More recently we have also introduced a wider range of Bi Polar Forceps and we have increased the range of standard off the shelf packs.

Since launching in 2004, our primary objective has been to supply the highest quality surgical instruments, in conjunction with first class customer service. Our high standards and attention to detail are reflected throughout the company, beginning with our manufacturing sites which adhere to ISO13485:2016, through to quality inspection, packaging, sterilisation, and the onward supply to the end user.

If you require any further information please call our customer service team on +44 121 386 8433. Alternatively to view our extensive product range visit us at www.blinkmedical.com.





blinkmedical



www.blinkmedical.com


Manufacturing Facility

All Blink Medical's surgical instruments are manufactured to the highest standard. Each of our manufacturing sites are ISO accredited and audited regularly to ensure compliance. All instruments undergo stringent inspection and cleaning protocols before being shipped to our state of the art distribution centre.



Logistics

Our UK Logistics Centre handles thousands of shipments every year, sending consignments to the UK and around the world. Whether ordering one box of forceps or a large consignment of customised sets, our highly trained professional staff are here to ensure your shipment is expertly handled and delivered on time.



Customer Service

Our customer service department will handle your order professionally and courteously, expediting it in a timely fashion. Should you have a query or require a quotation, our friendly, helpful staff will aim to deal with your request politely and quickly. Where possible we will process any quotation for a bespoke set the same day.

UK Quality Control

Our UK based quality control department is staffed by fully trained and experienced associates, who undertake the inspection of all instruments in accordance with the requirements of our ISO13485:2016 accredited Quality Management System.

All instruments are subject to stringent inspection in accordance with their intended use in a safe and effective manner. Only instruments meeting the acceptance criteria can be subsequently released for final packaging operations.

All instruments and procedure packs are assembled within a validated 'Class 8' cleanroom, with subsequent sterilisation of the packaged products completed ensuring sterility is achieved and maintained in accordance with the relevant regulatory requirements.


All production processes are subject to appropriate quality assurance control procedures ensuring that all product and regulatory requirements are met on a perpetual basis.



Internationally Available

Our exclusive International Distributor network covers the majority of Europe, Australia and New Zealand. Our main distributors carry the Ophthalmic range of instruments, however some also now carry the Cosmetic, ENT and Dermatology lines as well. Should you represent a company interested in becoming a Blink Distributor or alternatively if you are a potential customer wishing to locate your nearest Blink distributor, please visit our website at www.blinkmedical.com and click on the appropriate tab.

10-1415 - Forceps Birks Suture Tying
8mm Tying Platform, O/L 110mm




10-1302 - Forceps Utility Lister
1:2 Teeth 0.8mm with 6mm Tying Platform, O/L 114mm



10-1107 - Forceps Bonn
1:2 Teeth, Tying Platform 6mm, O/L 82mm


10-1106 - Forceps Bonn
1:2 Teeth, Tying Platform 6mm, O/L 105mm



10-9009 - Forceps Tying McPhearson
Straight Plain Tip 6mm Tying Platform O/L 85mm



11-7009 - Forceps Piers Collibri Notched
0.3mm Notch O/L 85mm



11-7010 - Forceps Collibri Toothed
1:2 Teeth 0.8mm Tooth O/L 65mm

11-7011 - Forceps Collibri Toothed
1:2 Teeth 0.8mm Tooth O/L 110mm



Ophthalmology

Tissue Forceps

10-1203 - Forceps St Martins

1:2 Teeth 0.6mm, O/L 90mm

10-1204 - Forceps St Martins

1:2 Teeth 0.3mm, O/L 90mm



11-6009 - Forceps Max Fine Curved

Smooth Tip, O/L 110mm



10-9008 - Forceps Tying McPhearson

Straight Plain Tip 6mm Tying Platform O/L 110mm



10-8008 - Forceps Hoskins Straight Notched

0.25mm Notched, 5mm Tying Platform, O/L 90mm

10-8009 - Forceps Hoskins Straight Notched

0.25mm Notched, 5mm Tying Platform, O/L 108mm



11-6016 - Forceps Tying Kelman McPhearson

5mm Tying Platform O/L 85mm

11-6015 - Forceps Tying Kelman McPhearson

5mm Tying Platform, O/L 105mm



11-6012 - Forceps Tying Curved

5mm Tying Platform, O/L 105mm



11-5051 - Forceps Capsulorhexis Utrata Curved

Flat Handle
Direct Action
Curved Shaft
0.4mm Sharp Tip
O/L 105mm



11-5050 - Forceps Capsulorhexis Utrata

Flat Handle
Direct Action
45 Degree Shaft
0.4mm Sharp Tip
O/L 105mm



11-5053 - Forceps Capsulorhexis Fulcher

Flat Handle
Direct Action
Curved Shaft
0.4mm Rounded Tip
O/L 105mm



11-5052 - Forceps Capsulorhexis

Flat Handle
Direct Action
45 Degree Shaft
O/L 105mm



11-5021 - Forceps Capsulorhexis Inamura

Flat Handle
Cross Action
Curved Shaft
0.4mm Sharp Tip
O/L 105mm



11-5011 - Forceps Capsulorhexis Inamura

Flat Handle
Cross Action
45 Degree Shaft
0.4mm Sharp Tip
O/L 105mm



Ophthalmology

Capsulorhexis Forceps

11-5013 - Forceps Capsulorhexis Utrata

Round Handle
Direct Action
Curved Shaft
0.4mm Sharp Tip
O/L 103mm



11-5012 - Forceps Capsulorhexis Utrata

Round Handle
Direct Action
45 Degree Shaft
0.4mm Sharp Tip
O/L 103mm



11-5023 - Forceps Capsulorhexis Micro Incision 1.8mm

Direct Action
Curved Shaft
0.4mm Sharp Tip
O/L 97mm




11-5027 - Forceps Capsulorhexis Micro Incision 1.8mm

Direct Action
45 Degree Shaft
0.4mm Sharp Tip
O/L 97mm



11-5020 - Forceps Capsulorhexis Inamura

Round Handle
Cross Action
Straight Shaft
0.4mm Sharp Tip
O/L 105mm




After listening to our customers and discussions with our quality, planning and production teams, we are now packaging our Capsulorhexis Forceps in rigid blister packs.




10-1300 - Forceps Epilation
Angled Tips
O/L 91mm




10-1301 - Forceps Suture Moorfields
16mm Serrated Tying Platform
O/L 107mm




10-1311 - Forceps Suture Moorfields Utility
14mm Serrated Tying Platform
O/L 107mm



11-1005 - Forceps Mosquito Curved
Curved Tip
O/L 100mm



10-1006 - Forceps Mosquito Straight
O/L 100mm



Ophthalmology

Scissors

20-7001 - Scissors Iris Straight
Rounded Tip O/L 100mm



20-7002 - Scissors Iris Straight
Sharp Tip O/L 100mm



20-6000 - Scissors Vannas Straight
O/L 95mm



21-6001 - Scissors Vannas Curved
O/L 95mm



20-4006 - Scissors Westcott Straight
Standard Tip O/L 120mm



20-4007 - Scissors Westcott Straight
Sharp Tip O/L 120mm



21-4008 - Scissors Westcott Curved
Sharp Tip O/L 120mm



21-4009 - Scissors Westcott Curved
Standard Tip O/L 120mm



21-4010 - Scissors Westcott Curved
Fine Tip O/L 120mm



21-5001 - Scissors Westcott Curved Small
Sharp Tip O/L 100mm



21-3010 - Scissors Ong
Sharp Tip O/L 92mm



21-4011 - Scissors Dolphin Nose
Dolphin Tip O/L 120mm




Do you perform Sub Tenon's Anaesthesia?

- > The NEW Conjunctival Probe™, a cost effective alternative to Westcott Scissors. The first dedicated instrument for making precise Sub Tenon incisions.
- > "The Conjunctival Probe™ would appear to be a safe and effective way of puncturing the conjunctiva to access the sub-Tenons space. It clearly appears to cause less conjunctiva trauma than Westcott scissors thereby minimising risk of conjunctival complications and reducing reflux of local anaesthetic"


Ophthalmic Anaesthesia Journal article June 2018

New




Introducing the New **SISTA**™ Technique Small Incision Sub Tenon Anaesthesia...


> The new **SISTA**™ technique using the revolutionary conjunctival probe for small incision sub tenon anaesthesia




1 Raise the Conjunctiva and Tenons using a pair of non-toothed Moorfields forceps, holding them at 90 degrees to the direction of the probe.




2 Holding the probe parallel to the surface of the sclera, puncture the Conjunctiva and Tenons.



3 Keeping the probe in contact with the sclera, follow the curve of the globe for approximately 1cm Sub Tenon.



4 Withdraw the probe but maintain the grip on the Conjunctiva and Tenons with the forceps.




5 Insert the Sub Tenon Cannula to inject the anaesthetic in the usual manner.

30-4000 - Speculum Barraquer
Solid Blade Angled Temporal L40mm W42mm



30-4001 - Speculum Adjustable
Sliding Angled Temporal L45mm



30-3001 - Speculum Burratto
Straight L57mm



30-3002 - Speculum Paediatric Spring
Solid Blade 7mm O/L35mm L30mm W30mm




31-3003 - Speculum Cook Adjustable
Solid Blades L40mm




31-3004 - Speculum Clarkes
Angled Temporal L80mm



31-3100 - Speculum Neo Natal
Angled Temporal W21mm O/L34mm



31-3101 - Speculum Neo Natal
Strong Spring Angled Temporal W21mm O/L34mm




31-3002 - Speculum Leibinger
Angled Temporal O/L65mm



Ophthalmology

Speculae


31-2003 - Speculum Barraquer
Angled Temporal L40mm W30mm



30-2004 - Speculum Barraquer Straight
L40mm W30mm




31-1006 - Speculum Kratz-Barraquer
Angled Temporal L40mm W30mm




30-1007 - Speculum Kratz-Barraquer
Straight L40mm W30mm





Temporal



Straight




31-3000 - Solid Blade Lebinger Speculum
Angled Temporal O/L65mm







40-5002 - Manipulator Lens Hook Sinskey Straight
Tip L 0.5mm O/L 124mm




41-4003 - Manipulator Phaco Chopper Mackool
Tip Angled at 45 Degrees
Tip L 2mm O/L 120mm




41-5001 - Manipulator Lens Hook Sinskey Angled
Tip Angled at 45 Degrees
Tip L 0.5mm O/L 120mm



41-7008 - Chopper Greens Phaco
Tip Angled at 45 Degrees
Tip L 10mm O/L 117mm




41-3004 - Manipulator Spatula
Tip Angled at 45 Degrees
Tip L 12mm W0.75mm O/L 135mm




41-7009 - Chopper Claw
Tip Angled at 45 Degrees
Tip L 10mm O/L 123mm



41-3005 - Neuhann Spatula
Tip Angled at 45 Degrees
Tip Curved Angled



41-6007 - Paddle Manipulator Drysdale
Tip Angled at 45 Degrees
Tip L 11mm Spatula 1mm O/L 120mm



41-2005 - Rotator Nucleus Y Shaped
On horizontal axis
Tip Angled at 45 Degrees
Tip L 10mm O/L 124mm



41-4005 - Chopper Sadiq
Tip Angled at 45 Degrees
Tip L 11mm O/L 121mm




41-1006 - Lester IOL Manipulator
Tip Angled at 45 Degrees
Tip L 12mm 1mm Branch O/L 125mm



41-4007 - Haefliger Chopper
Tip Angled at 45 Degrees
Tip L 11mm O/L 115mm



41-8009 - Manipulator Mushroom Tipped
Tip Angled at 45 Degrees
Tip L 11mm O/L 123mm





41-4006 - Spatula Govan
Tip Angled at 45 Degrees
Tip L 10mm O/L 130mm



Ophthalmology



41-9010 - Manipulator Sinskey/Mushroom Double Ended

Tips Angled at 45 Degrees
Sinskey Tip Length 0.65mm, Mushroom Tip
O/L 123mm




41-9011 - Manipulator Aycliffe/Sadiq Double Ended Chopper

Tips Angled at 45 Degrees
Aycliffe Tip 0.5mm and Sadiq Tip 1.5mm




41-9012 - Manipulator Mushroom/Sadiq Double Ended

Tips Angled at 45 Degrees
Sadiq Tip, Mushroom Tip 1.5mm



Refractive

51-0000 - Flap Lifter Lasik - Tip L12mm O/L 122mm



51-0001 - Epithelial Spatula - Spatula L7mm O/L 113mm



51-0002 - Alcohol Well - I/D 9MM O/L 114mm



90-2008 - Probe Liebricht 000+00
O/L125mm



90-2009 - Probe Liebricht 00+0
O/L125mm



90-2010 - Probe Liebricht 1+2
O/L125mm



90-2011 - Probe Liebricht 3+4
O/L125mm



90-1011 - Dilator Punctal Double Ended
O/L138mm



90-1012 - Nettleships Dilator
Tip 20mm O/L118mm



90-8012 - Bulldog Clamp
O/L 52mm



90-9013 - Desmares Retractor
14mm Width O/L130mm



Ophthalmology

Oculoplastic

90-7001 - Needle Holder Castroviejo
Straight with lock 12mm Tying Platform O/L135mm



90-7002 - Needle Holder Castroviejo
Without lock 12mm Tying Platform O/L135mm




91-7003 - Needle Holder Castroviejo
Delicate Curve 12mm Tying Platform O/L135mm



91-7004 - Needle Holder Troutman
6mm Tying Platform O/L105mm



90-5004 - Hook Stabismus Graefe Medium
Tip L10mm O/L140mm



90-5005 - Hook Alexandra Muscle
Tip L11mm O/L130mm



90-4005 - Clamp Chalazion Small
Clamp Diameter 11mm O/L90mm



90-4006 - Clamp Chalazion Medium
Clamp Diameter 18mm x11mm O/L80mm



90-4007 - Clamp Chalazion Large
Clamp Diameter 23mm x15mm O/L90mm



90-4008 - Desmares Entropian Clamp Left
Clamp Diameter 24mm x 16mm O/L95mm




90-4009 - Desmares Entropian Clamp Right
Clamp Diameter 24mm x 16mm O/L95mm




90-3007 - Curette Double Ended
2mm & 3mm cup O/L130mm



90-3008 - Curette Double Ended
2mm & 5mm cup O/L130mm



90-3009 - Curette Double Ended
2mm & 4mm cup O/L130mm




Ophthalmology

For Intravitreoretinal Injections

90-1014 - Caliper and Rule



90-1015 - Double Ended Marker
3.5 - 4.0mm O/L 113mm



90-1016 - Marker Fixator
3.5mm O/L 70mm



90-1017 - Braunstien Marker
3.5 - 4.0mm O/L 80mm



Other

90-1020 - Thornton Ring



91-6001 - Alignment Marker, 3, 6, 9 O'clock.



91-6000 - Degree Gauge 11.7mm I.D.



91-6002 - Toric IOL Marker 11mm



Ophthalmology

Standard Ophthalmic Sets

001 - Sub Tenon Set - Standard Tip 001-A - Sub Tenon Set - Fine Tip



Qty	Code	Description
001		
1	21-4009	Scissors Westcott Curved Standard Tip
1	31-1006	Speculum Kratz-Barraquer Angled Temporal
1	10-1301	Forceps Suture Moorfields
1		Wrap
001A		
1	21-4010	Scissors Westcott Curved Fine Tip
1	31-1006	Speculum Kratz-Barraquer Angled Temporal
1	10-1301	Forceps Suture Moorfields
1		Wrap

270 - CP Sub Tenon Basic Set



Qty	Code	Description
1	10-1301	Forceps Suture Moorfields
1	31-1006	Speculum Kratz-Barraquer Angled Temporal
1	CP01	Conjunctival Probe
1		Wrap

300 - Complete Sub Tenon Set




Qty	Code	Description
1	10-1301	Forceps Suture Moorfields
1	31-2003	Speculum Barraquer Angled Temporal
1	21-4010	Scissors Westcott Curved Fine Tip
1		Sub Tenon Cannula
1		Syringe 5ml Luer Lock
5		Swab Non Woven 10x10cm
1		Tray with Gallipot Holder
1		Wrap

290 - Complete CP Sub Tenon Set



Qty	Code	Description
1	10-1301	Forceps Suture Moorfields
1	31-2003	Speculum Barraquer Angled Temporal
1	CP01	Conjunctival Probe
1		Sub Tenon Cannula
1		Syringe 5ml Luer Lock
5		Swab Non Woven 10x10cm
1		Tray with Gallipot Holder
1		Wrap

023 - Basic Intravitreal Injection Set



Marker double ended 3.5mm-4mm

Qty	Code	Description
1	90-1015	Marker double ended 3.5mm-4mm
1	10-1104	Forceps Suture Bonn 1:2 Teeth
1	30-1007	Speculum Kratz-Barraquer Straight
1		Wrap

208-A - Complete Ophthalmic Intravitreal Set



Qty	Code	Description
1	31-2003	Speculum Barraquer Angled Temporal
1	90-1015	Marker Double Ended Scleral
3		Eye Spears
3		Double ended Cotton Buds
1		Plastic Rampley Sponge Holders
1		Drape 40x40cm Slit Incise Drape
1		60ml Gallipot
1		Swab Non Woven 10x10cm
1		Tray with Gallipot Holder
1		Wrap

092 - Chalazion Set - Small Clamp

111 - Chalazion Set - Medium Clamp



Qty	Code	Description
092		
1	90-4005	Chalzion Clamp Small
1	90-3009	Curette Double Ended Small/Medium
1		No11 Scalpel
1		Wrap
111		
1	90-4006	Chalzion Clamp Medium
1	90-3009	Curette Double Ended Small/Medium
1		No11 Scalpel
1		Wrap

210 - Ophthalmic Suture Removal Set



Qty	Code	Description
1	10-9007	Forceps Tying McPhearson Straight
1	21-6004	Scissors Vannas Curved
1		Wrap

Custom Packs

Efficiency in how surgical procedures are prepared for and carried out is of ever increasing importance to our customers...

Hospitals are increasingly being focussed on efficiency savings within all areas, with a high focus on their Theatre Departments. To support this drive for greater efficiency, Blink Medical have developed a range of standard and Custom Procedure Packs.

Our sales team are trained to identify customer needs and suggest which option of Standard or Custom Packs best meets their requirements. The scope of our Custom Procedure Pack offering covers surgical areas such as Ophthalmology, Dermatology, ENT and Cosmetic Surgery.



- Increased Efficiency
- Improved Infection Control
- Cost Control
- Improved Product Traceability
- Clinical Waste Reduction



Page 1 - Introduction

Page 3 - Company Information

Page 6 - Ophthalmology Tissue Forceps

10-1415 - Forceps Birks Suture Tying
 10-1302 - Forceps Utility Lister
 10-1107 - Forceps Bonn
 10-1106 - Forceps Bonn Micro
 10-9009 - Forceps Tying McPhearson
 11-7009 - Forceps Pierse Colibri Notched
 11-7010 - Forceps Colibri Toothed
 11-7011 - Forceps Colibri Toothed

Page 7 - Ophthalmology Tissue Forceps

10-1203 - Forceps St Martins
 10-1204 - Forceps St Martins
 11-6009 - Forceps Max
 10-9008 - Forceps Tying McPhearson
 10-8008 - Forceps Hoskins Straight Notched
 10-8009 - Forceps Hoskins Micro
 11-6016 - Forceps Tying Kelman McPhearson
 11-6015 - Forceps Tying Kelman McPhearson
 11-6012 - Forceps Tying Curved

Page 8 - Ophthalmology Capsulorhexis Forceps

11-5011 - Forceps Capsulorhexis Inamura
 11-5021 - Forceps Capsulorhexis Inamura
 11-5050 - Forceps Capsulorhexis Ultrata
 11-5051 - Forceps Capsulorhexis Ultrata Curved
 11-5052 - Forceps Capsulorhexis
 11-5053 - Forceps Capsulorhexis Fulcher

Page 9 - Ophthalmology Capsulorhexis Forceps

11-5013 - Forceps Capsulorhexis Ultrata
 11-5012 - Forceps Capsulorhexis Ultrata
 11-5023 - Forceps Capsulorhexis Micro Incision
 11-5027 - Forceps Capsulorhexis Micro Incision

Page 10 - Ophthalmology Forceps Other

10-1300 - Forceps Epilation
 10-1301 - Forceps Suture Moorfields
 10-1311 - Forceps Suture Moorfields Utility
 11-1005 - Forceps Mosquito Curved
 11-1006 - Forceps Mosquito Straight

Page 11 - Ophthalmology Scissors

20-7001 - Scissors Iris Straight
 20-7002 - Scissors Iris Straight
 20-6000 - Scissors Vannas Straight
 21-6001 - Scissors Vannas Curved
 20-4006 - Scissors Westcott Straight Blunt
 20-4007 - Scissors Westcott Straight Sharp

Page 12 - Ophthalmology Scissors

21-4008 - Scissors Westcott Curved Sharp
 21-4009 - Scissors Westcott Curved Blunt
 21-4010 - Scissors Westcott Curved Blunt
 21-5001 - Scissors Westcott Curved Small
 21-3010 - Scissors Ong
 21-4011 - Scissors Dolphin Nose

Page 13 - Conjunctival Probe

Page 14 - Ophthalmology Speculae

30-4000 - Speculum Barraquer
 30-4001 - Speculum Adjustable
 30-3001 - Speculum Burratto
 30-3002 - Speculum Paediatric Spring
 31-3003 - Speculum Cook Adjustable
 31-3004 - Speculum Clarkes
 31-3100 - Speculum Neo Natal
 31-3101 - Speculum Neo Natal
 31-3002 - Speculum Leibinger

Page 15 - Ophthalmology Speculae

31-2003 - Speculum Barraquer
 30-2004 - Speculum Barraquer Straight
 31-1006 - Speculum Kratz-Barraquer
 30-1007 - Speculum Kratz-Barraquer
 31-3000 - Speculum Leibinger

Page 16 - Ophthalmology Lens Manipulation

40-5002 - Manipulator Lens Hook Sinskey Straight
 41-5001 - Manipulator Lens Hook Sinskey Angled
 41-3004 - Manipulator Spatula
 41-2005 - Rotator Nucleus Y Shaped
 41-1006 - Lester IOL Manipulator
 41-3005 - Neuharn Spatula
 41-3009 - Manipulator Mushroom Tipped
 41-4006 - Spatula Govan
 41-4003 - Manipulator Phaco Chopper Mackool
 41-7008 - Chopper Greens Phaco
 41-7009 - Chopper Claw
 41-6007 - Padille Manipulator Drysdale
 41-4005 - Chopper Sadiq
 41-4007 - Haeltiger Chopper

Page 17 - Ophthalmology Lens Manipulation

41-9010 - Manipulator Sinskey/Mushroom
 41-9011 - Manipulator Aycliffe/Chopper
 41-9012 - Manipulator Mushroom/Sadiq
 51-0000 - Flap Lifter Lasik
 51-0001 - Epithelial Spatula
 51-0002 - Alcohol Well

Page 18 - Ophthalmology Oculoplastic

90-2008 - Probe Liebricht 000+00
 90-2009 - Probe Liebricht 00+0
 90-2010 - Probe Liebricht 1+2
 90-2011 - Probe Liebricht 3+4
 90-1011 - Dilator Punctal
 90-1012 - Nettleships Dilator
 90-8012 - Bulldog Clamp
 90-9013 - Desmares Retractor

Page 19 - Ophthalmology Oculoplastic

90-7001 - Needle Holder Castroviejo
 90-7002 - Needle Holder Castroviejo
 91-7003 - Needle Holder Castroviejo
 90-7004 - Needle Holder Troutman
 90-5004 - Hook Stabismus Graefe Medium
 90-5005 - Hook Alexandra Muscle

Page 20 - Ophthalmology Oculoplastic

90-4005 - Clamp Chalazion Small
 90-4006 - Clamp Chalazion Medium
 90-4007 - Clamp Chalazion Large
 90-4008 - Desmares Entropion Clamp Left
 90-4009 - Desmares Entropion Clamp Right
 90-3007 - Curette Double Ended
 90-3008 - Curette Double Ended
 90-3009 - Curette Double Ended

Page 22 - Ophthalmology IntraVitreoretinal Injections

90-1014 - Caliper and Rule
 90-1015 - Double Ended Marker
 90-1016 - Marker Fixator
 90-1017 - Braunsten Marker

Page 22 - Ophthalmology Other

91-6000 - Degree Gauge 11,7mm I.D.
 91-6001 - Alignment Marker, 3, 6, 9 O'clock
 90-1020 - Thornton Ring
 91-6002 - Toric IOL Marker 11mm

Page 23 - Ophthalmology Standard Ophthalmic Sets

001 - Sub Tenon Set
 001A - Sub Tenon Set
 270 - CP Sub Tenon Basic Set
 300 - Complete Sub Tenon Set
 290 - Complete CP Sub Tenon Set

Page 24 - Ophthalmology Standard Ophthalmic Sets

023 - Basic IntraVitreoretinal Injection Set
 208A - Complete Ophthalmic IntraVitreoretinal Set
 092 - Chalazion Set
 111 - Chalazion Set
 210 - Ophthalmic Suture Removal Set

Page 25 - Custom Packs

Page 26 - General Surgery Standard Sets

095 - Punch Biopsy Set
 185 - Biopsy Full Set
 074 - Excision Set
 039 - Excision Set Full

Page 27 - General Surgery Standard Sets

038 - Minor Dermatology Set (A)
 020 - Minor Ops Set (B)
 047 - Minor Ops Set Full
 263 - Carpel Tunnel Set

Page 28 - General Surgery Standard Sets

066 - Suture Removal Set
 252 - Perio Probe Dental Set

A KATENA BRAND

blink^omedical



twitter.com/BlinkMedical



www.facebook.com/blinkmedical



www.linkedin.com/company/blink-medical/

Blink Medical Ltd

Radway Road
Shirley
Solihull
B90 4NS

Tel: **+44 (0) 121 386 8433**

Fax: **+44 (0) 121 277 4949**

Web: **www.blinkmedical.com**



Version 5.0

DISTRIBUTORE ITALIANO

SIR Oftalmica Srl

Via Torriani 19/C

22100 - COMO

Tel. : 031-570869

info@siroftalmica.com

www.siroftalmica.com

P.IVA : 01814920136

Cell. A. Molteni: 328.1808176

Cell. L. Ronchi: 335.6167324



SIRoftalmica