



## Suggested Basic Lid Set




This set provides the surgeon with an array of instruments that will allow for a variety of eyelid surgical procedures including the correction of ectropia and entropia, as well as ptosis repairs and eyelid trauma. The set also includes a curette and clamp for removing chalazia.

### 6-664 Otaka Lid Plate

- 23.5mm wide
- 20.4mm wide
- Length 55mm




### 2-401NR8 Straight 0.65mm Notched Round Handle Forceps




- Pierce 0.65mm notched, 6mm tying platforms
- Straight shafts
- 8mm diameter round handle, length 115mm

### 2-116N Straight Toothed Forceps (0.3mm)




- 0.3mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- New long flat handle, length 115mm

### 1-500 Westcott Style Stitch Scissors (straight)



- Semi-sharp tips, straight blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 116mm


### 1-512 Westcott Style Tenotomy Scissors (Right Handed)



- Sharp tips, curved blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 116mm




### 6-805 Bulldog Clip



- 8.5mm serrated cross action jaws
- Length 46mm




### 6-628 Desmarres Lid Retractor (size 0)



- Blade width 11mm
- Flat handle, length 134mm




### 3-202 DK Needle Holder (curved with lock)



- 9mm medium curved jaws
- Round handle with lock, tag spring
- Length 110mm



### 3-203 DK Needle Holder (curved)




- 9mm medium curved jaws
- Round handle tag spring
- Length 108mm




### 2-640 DK Desmarres Chalazion Forceps / Clamp

- 11.5mm x 17.5mm inside oval circular working diameter upper plate
- 14.0mm x 20mm diameter lower plate
- Polished finish lower plate
- Standard handle, locking thumb screw, length 88mm




### 6-641 DK Meyerhoefer Chalazion Curette, 1mm cup




- 1mm cup
- Straight shaft
- Round Handle, length 125.5mm

### T7010 Sterilising Tray (for 10 instruments)

- External dimensions: 26.4cm x 16.2cm x 2.6cm
- Internal dimensions: 25.4cm x 15.2cm x 1.8cm (without silicone mat)
- Suitable for 10 instruments
- One silicone mat



Visit : [www.duckworth-and-kent.com](http://www.duckworth-and-kent.com) to view  
a wide range of ophthalmic titanium instruments



As well as product information for a wide range of instruments, you will also find detailed literature and informative videos showing some of our instruments in action.

## Duckworth & Kent Promotional Ranges

*See the Quality Feel the Difference*

**CORE VALUES**  
*Collection*

**E-RANGE**  
*Instruments*

*Essential*  
**CATARACT SET**

Connect



[www.duckworth-and-kent.com](http://www.duckworth-and-kent.com)

7 Marquis Business Centre  
Royston Road, Baldock  
Herts SG7 6XL England

Tel: +44 (0)1462 893254  
Fax: +44 (0)1462 896288  
Email: [info@duckworth-and-kent.co.uk](mailto:info@duckworth-and-kent.co.uk)



© November 2014 Duckworth & Kent  
04.07.2016

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.

*at the Leading Edge*