




Suggested Consulting Room - Eye Emergency Set




An invaluable set of instruments which assists the surgeon in performing procedures which may be carried out in the emergency eye room or consulting room, such as epilation, removal of corneal and eyelid sutures, removal of contact lenses and basic repair of eyelid and conjunctival traumatic wounds.

9-556 *Kratz Barraquer Speculum*

- 14.5mm open blades
- Angled to rest temporally
- Single piece construction




2-510-1E *Moorfields Utility Forceps*




- Delicate serrated tips, 1mm wide
- Straight shafts
- Flat handle, length 115mm
- Design registration number 004383396-0001




1-501B *DK Westcott Style Stitch Scissors (curved)*



- Blunt tips, curved blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 115mm




2-110E *Bonn Suturing Forceps (0.12mm)*



- 0.12mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- Flat handle, length 90mm
- Design registration number 004383396-0001




1-118 *DK Castroviejo Miniature Corneal Scissors (curved)*



- Blunt tips, curved blades
- Universal cutting blades
- Cut length 6mm, tip to pivot length 11mm
- Standard handle, length 116mm




2-501N Curved Tying Forceps



- 6mm tying platforms
- Curved shafts
- Standard handle, length 83mm




3-205 Fine Curved Needle Holder



- 7mm delicate curved jaws
- Round handle, tag spring
- Length 107mm



2-2195E Cilia Forceps




- Tapered, smooth tips
- 30° angled shafts, tip to angle length 7mm
- Short ridged handle, length 90mm
- Design registration number 004383396-0001



T7010 Sterilising Tray (suitable for 10 instruments)

- External dimensions: 26.4cm x 16.2cm x 2.6cm
- Internal dimensions: 25.4cm x 15.2cm x 1.8cm (without silicone mat)
- Suitable for 10 instruments
- One silicone mat



www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.co.uk



© March 2011 Duckworth & Kent
20.02.18

DAK® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.

at the Leading Edge