




Epilation, Meibomian & Chalazion Set




A comprehensive set of instruments designed to allow the surgeon to perform epilation, expression of meibomian glands and safe and effective incision and curettage of eyelid chalazia. Our full range of curette sizes is provided.

2-635 Kudo Meibomian Gland Pressing Forceps



- 60° angled shafts, tip to angle length 11.5mm
- Internal faces of tip are cross hatched
- Standard cross action handle, length 112mm




2-640 DK Desmarres Chalazion Forceps / Clamp

- 11.5mm x 17.5mm inside oval circular working diameter upper plate
- 14.0mm x 20.0mm diameter lower plate
- Polished finish lower plate
- Standard handle, locking thumb screw, length 88.0mm




6-641 DK Meyerhoefer Chalazion Curette, 1mm cup



- 1mm cup
- Straight shaft
- Round Handle, length 125.5mm




6-641-1 DK Meyerhoefer Chalazion Curette, 1.5mm cup



- 1.5mm cup
- Straight shaft
- Round Handle, length 125.5mm




6-641-2 DK Meyerhoefer Chalazion Curette, 2mm cup



- 2mm cup
- Straight shaft
- Round Handle, length 125.5mm




6-641-3 DK Meyerhoefer Chalazion Curette, 2.5mm cup




- 2.5mm cup
- Straight shaft
- Round Handle, length 125.5mm

6-641-4 DK Meyerhoefer Chalazion Curette, 3mm cup



- 3mm cup
- Straight shaft
- Round Handle, length 125.5mm


2-2195E Cilia Forceps



- Tapered, smooth tips
- 30° angled shafts, tip to angle length 7mm
- Short ridged handle, length 90mm
- Design registration number 004383396-0001

T7010 Sterilising Tray (for 10 instruments)

- External dimensions: 26.4cm x 16.2cm x 2.6cm
- Internal dimensions: 25.4cm x 15.2cm x 1.8cm (without silicone mat)
- Suitable for 10 instruments
- One silicone mat



www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.co.uk



© May 2014 Duckworth & Kent
20.02.18

DAK® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.

at the Leading Edge