




## Suggested Lacrimal Set



A collection of instruments which provide the surgeon with a comprehensive range of dilators for enlarging the lacrimal punctum and performing 3-snip punctoplasties, as well as our complete array of lacrimal probes for nasolacrimal drainage investigation and treatment.


### 6-180 Lacrimal Dilator



- 0.1mm diameter pointed tip
- 10° inclusive angle
- Round handle, length 122mm




### 6-180-1 Lacrimal Dilator



- 0.2mm diameter rounded tip
- 10° Inclusive Angle
- Round Handle, length 122mm




### 6-181 Lacrimal Dilator




- 0.45mm diameter blunt tip
- Round Handle, length 120mm



### 6-182-2 Otaka Dilator




- Bullet shape tip to ease insertion
- Double ended instrument
- Marks at diameters 0.5, 0.6 and 0.7mm



- Marks at diameters 0.8, 0.9, 1.0, 1.1 and 1.2mm
- Round handle, length 103mm


### 2-100E Pierse Notched Forceps



- Pierce 0.25mm notched, 6mm tying platforms
- Straight shafts
- Flat handle, length 90mm
- Design registration number 004383396-0001




### 1-110 Vannas Scissors (straight)



- Sharp tips, straight blades
- Cut length 5mm, tip to pivot length 9mm
- Standard handle, length 95mm




### 6-656 DK Lacrimal Probe




- High quality surface finish over first 5mm of tips
- Overall length 130mm
- Tip sizes: 0000 and 000




### 6-656-1 DK Lacrimal Probe




- High quality surface finish over first 5mm of tips
- Overall length 130mm
- Tip sizes: 00 and 0




### 6-656-2 DK Lacrimal Probe



- High quality surface finish over first 5mm of tips
- Overall length 130mm
- Tip sizes: 1 and 2




### 6-656-3 DK Lacrimal Probe





- High quality surface finish over first 5mm of tips
- Overall length 130mm
- Tip sizes: 3 and 4

### 6-656-4 DK Lacrimal Probe



- High quality surface finish over first 5mm of tips
- Overall length 130mm
- Tip sizes: 5 and 6


### 6-656-5 DK Lacrimal Probe




- High quality surface finish over first 5mm of tips
- Overall length 130mm
- Tip sizes: 7 and 8

### T7010 Sterilising Tray (for 10 instruments)

- External dimensions: 264mm x 162mm x 25mm
- Suitable for 10 instruments
- One silicone mat



Visit : [www.duckworth-and-kent.com](http://www.duckworth-and-kent.com) to view  
over 800 ophthalmic titanium instruments



As well as product information for over 800 instruments, you will also find detailed literature and informative videos showing some of our instruments in action.

## Duckworth & Kent Promotional Ranges

*See the Quality Feel the Difference*

**CORE VALUES**  
*Collection*

**E-RANGE**  
*Instruments*

*Essential*  
**CATARACT SET**

Connect



[www.duckworth-and-kent.com](http://www.duckworth-and-kent.com)

7 Marquis Business Centre  
Royston Road, Baldock  
Herts SG7 6XL England

Tel: +44 (0)1462 893254  
Fax: +44 (0)1462 896288  
Email: [info@duckworth-and-kent.co.uk](mailto:info@duckworth-and-kent.co.uk)



© March 2011 Duckworth & Kent  
05.12.17

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.

*at the Leading Edge*