



## Suggested LASEK Set



A set providing the surgeon with instruments to assist in LASEK procedures, including a speculum as well trephines and a variety of epithelial edge lifters.

### 9-588-3 Lieberman Adjustable Speculum


- 15mm open blades
- Angled to rest temporally
- Adjustable with thumb screw



### 6-946 DK 9.0mm Trephine Guide - Alcohol Chamber



VIEW ON ARROW




- Alcohol chamber, height 6mm
- Internal diameter 0.5mm larger than incision
- Designed as a guide with the Ø9mm incision DK Epithelial Trephine 6-926
- Round handle, length 125mm

### 6-947 DK 9.5mm Trephine Guide - Alcohol Chamber



VIEW ON ARROW



- Alcohol chamber, height 6mm
- Internal diameter 0.5mm larger than incision
- Designed as a guide with the Ø9.5mm incision DK Epithelial Trephine 6-927
- Round handle, length 125mm

### 6-926 DK Epithelial Trephine for LASEK




- 9mm
- Used with 6-946, DK 9mm trephine guide - alcohol

## 6-866 DK Epithelial Separator / Lifter




- 5mm triangular-shaped tip
- 45° angled shaft
- Round handle, length 120mm

## 6-103-1 DK Double Ended Spatula (for repositioning epithelial flap)




- 0.8mm diameter, double ended
- 45° angled shaft, tip to angle length 12mm, curved shaft
- Round handle, length 127mm



## T7010 Sterilising Tray (for 10 instruments)

- External dimensions: 26.4cm x 16.2cm x 2.6cm
- Internal dimensions: 25.4cm x 15.2cm x 1.8cm (without silicone mat)
- Suitable for 10 instruments
- One silicone mat



[www.duckworth-and-kent.com](http://www.duckworth-and-kent.com)

7 Marquis Business Centre    Tel: +44 (0)1462 893254  
Royston Road, Baldock        Fax: +44 (0)1462 896288  
Herts SG7 6XL England        Email: [info@duckworth-and-kent.co.uk](mailto:info@duckworth-and-kent.co.uk)



© March 2011 Duckworth & Kent  
16.03.18

DAK® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.

*at the Leading Edge*