




## Retinal Detachment Set



A set of instruments compiled to allow the surgeon to perform conventional extra-ocular retinal detachment surgery using plumbs, buckles and encircling bands to secure the retina in place.


### 9-551 Barraquer Adult Speculum

- 14.5mm closed blades
- Angled to rest temporally
- Single piece construction




### 9-573 Barraquer Paediatric Speculum

- 9mm closed blades
- Angled to rest temporally
- Single piece construction




### 2-108 Rabkin Blepharoplasty Tissue Forceps (0.5mm)



- 0.5mm, 1 x 2 flat teeth, 6mm tying platforms
- Straight shafts
- Flat handle, length 116mm




### 2-500-4E Chihara Conjunctival Forceps



- 6mm platforms
- Serrated tips
- Straight shafts
- Flat handle, length 90mm
- Design registration number 004383396-0001




### 6-805 Bulldog Clip



- 8.5mm serrated cross action jaws
- Length 46mm




### 1-501 Westcott Style Stitch Scissors (curved)



- Semi-sharp tips, curved blades
- Cut length 10mm, tip to pivot length 20mm
- Flat handle, length 116mm




### 1-500 Westcott Style Stitch Scissors (straight)



- Semi-sharp tips, straight blades
- Cut length 10mm, tip to blade length 20mm
- Flat handle, length 115mm




### 1-500B DK Westcott Style Stitch Scissors (straight)



- Blunt tips, straight blades
- Cut length 9mm, tip to pivot length 20mm
- Flat handle, length 115mm





### 6-625 Jameson Muscle Hook



- 2mm bulbous tip, tip to curve length 8mm
- Round handle, length 136mm




### 6-635-3 Koura Double Ended Scleral Depressor



- 2mm round and 2mm cylinder tips
- 2x slits down side, 0.5mm diameter x 6mm length
- Round handle, length 124mm

The slits down the side of the tip aid in gripping the sclera when depressed.


### 2-510-1E Moorfields Utility Forceps




- Delicate serrated tips, 1mm wide
- Straight shafts
- Flat handle, length 115mm
- Design registration number 004383396-0001

### T7020 Sterilising Tray (suitable for 20 instruments)

- External dimensions: 395mm x 266mm x 25mm
- Suitable for 20 instruments
- One silicone mat



Visit : [www.duckworth-and-kent.com](http://www.duckworth-and-kent.com) to view  
over 800 ophthalmic titanium instruments



As well as product information for over 800 instruments, you will also find detailed literature and informative videos showing some of our instruments in action.

## Duckworth & Kent Promotional Ranges

*See the Quality Feel the Difference*

**CORE VALUES**  
*Collection*

**E-RANGE**  
*Instruments*

*Essential*  
**CATARACT SET**

Connect



[www.duckworth-and-kent.com](http://www.duckworth-and-kent.com)

7 Marquis Business Centre  
Royston Road, Baldock  
Herts SG7 6XL England

Tel: +44 (0)1462 893254  
Fax: +44 (0)1462 896288  
Email: [info@duckworth-and-kent.co.uk](mailto:info@duckworth-and-kent.co.uk)



© August 2017 Duckworth & Kent  
12.12.17

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.

*at the Leading Edge*