




Suggested Veterinary Set




A set that provides the veterinary eye surgeon with instruments to assist in intra-ocular surgery.

9-559 DK Speculum

- 14mm large closed blades
- Curved to rest temporarily
- Single piece construction




2-110E Bonn Suturing Forceps (0.12mm)



- 0.12mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- Flat handle, length 90mm
- Design registration number 004383396-0001




2-132N Colibri Toothed Forceps (0.12mm)



- 0.12mm, 1 x 2 teeth, tip length 2mm
- 6mm tying platforms
- Colibri style shafts
- Flat handle, length 83mm



4-115 DK Retractable Diamond Knife, 1.0mm Trifacet




W=1 T=0.17 L=3.5

- Retractable long handle
- Length 123.5mm




9-650 DK Castroviejo Style Marking Caliper



- 0.25mm x 0.25mm delicate marking tips
- Marks 0mm to 15mm in 0.25 increments
- Marks from centre of tips
- Adjustable thumb screw
- Standard caliper handle




1-400 *Castroviejo Corneal Scissors (curved to right)*



- Blunt tips, curved blades, right
- Cut length 5mm, tip to pivot length 11mm
- Standard handle, length 95mm




1-401 *Castroviejo Corneal Scissors (curved to left)*



- Blunt tips, curved blades, left
- Cut length 5mm, tip to pivot length 11mm
- Standard handle, length 95mm



2-716GER8-1 *Inamura Round Handle Capsulorhexis Forceps*




- Sharp pointed serrated tips
- Cross action tips, 1.5mm width at pivot box
- Curved shaft, tip to pivot point 10.5mm
- Marks on shaft 2.5mm and 5mm denote desired size of capsulorhexis



- Tips angled 45° from shaft
- Round handle, overall length 111mm
- Design registration number 004383396-0001

2-2-819 *Masaoka Paddle Prechopper Forceps*




Tip height 1mm

- Tips fully open to 2.5mm
- Arrow tip shape for penetrating nucleus
- Flattened paddle area of tips for nucleus splitting, tip height 1mm



- 1.8mm maximum width at incision point
- Straight shafts, reverse cross action style
- Round handle, length 123mm



8-602 *Mackool® Hydrodissection Cannula*



- 23 gauge tube (Ø0.65mm), 0.9mm x 0.3mm flattened tip
- 30° angled shaft, tip to angle length 8mm




6-083-1 Hara Nucleus Divider (curved)



- Curved, sharp inner sides of tip, tip length 1.1mm
- Cutting edge 35° to axis
- 45° angled shaft, tip to angle length 11mm
- Round handle, length 115mm




8-635-3 Ogawa I/A Cannula



- 18 gauge, thin wall, 0.8mm irrigation port, 0.3mm aspiration port
- 30° angled shaft, tip to angle length 7mm
- Coaxial irrigation and aspiration
- Irrigation through main hub and aspiration through side port
- Silicone tubing and luer fitting supplied




6-099 Mackool Barraquer Spatula



- 0.45mm x 0.25mm
- Angled shaft, tip to angle length 15mm
- Round handle, length 116mm




6-250 DK Sinskey Hook



- 0.18mm diameter tip pointing down
- Angled shaft, tip to angle length 10mm
- Round handle, length 119mm




6-500 Fenzl Hook



- 0.13mm dialler
- Angled shaft, tip to angle length 9mm
- Round handle, length 117mm




2-500E McPherson Tying Forceps




- 6mm tying platforms
- Straight shafts
- Flat handle, length 90mm
- Design registration number 004383396-0001

3-206 Nordan Needle Holder (curved with lock)




- 7mm delicate curved jaws
- Round handle with lock, tag spring
- Length 107mm

33-220 DK Needle Holder (Straight with lock)



- 9mm medium straight jaws, blunt tip
- Round handle with lock, tag spring
- Length 136mm


1-121 DK Vannas Scissors (curved)



- Sharp tips, curved blades
- Cut length 9mm, tip to pivot length 15mm
- Standard handle, length 101mm

T7010 Sterilising Tray (for 10 instruments)

- External dimensions: 26.4cm x 16.2cm x 2.6cm
- Internal dimensions: 25.4cm x 15.2cm x 1.8cm (without silicone mat)
- Suitable for 10 instruments
- One silicone mat



www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.co.uk

© March 2011 Duckworth & Kent
05.12.17



D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeon's names or by the category of the instrument.

at the Leading Edge