




LRI Set





A set of instruments that allows the surgeon to perform limbal relaxing incisions and thereby reduce or eliminate corneal astigmatism.

9-556 *Kratz Barraquer Speculum*

- 14.5mm open blades
- Angled to rest temporally
- Single piece construction




9-510R *DK Fine Thornton Fixation Ring*



- 3/4 ring, 13mm diameter ring with nine point fixation
- Pivot swivel round handle, length 112mm


9-729 *DK Axis Marker (2 blades)*



- 2 axial blades
- 4mm inside diameter, 11.6mm outside diameter
- Lowest profile with centre pointer
- Round handle, length 124mm

— —
Marks created

9-706 *Packard-Rosen LRI Degree Marker / Fixation*




VIEW ON ARROW



- Marks 0° to 90° in 10° increments x 4
- 12mm internal ring diameter
- 12 point fixation
- Flat handle, length 117mm


2-110E Bonn Suturing Forceps (0.12mm)



- 0.12mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- Flat handle, length 90mm
- Design registration number 004383396-0001



4-620 Wallace LRI Diamond Knife (600 micron preset blade depth)




- 0.2mm front flat tip
- Blade preset at 600 micron
- Retractable handle, length 120mm



T7010 Sterilising Tray (suitable for 10 instruments)

- External dimensions: 264mm x 162mm x 25mm
- Suitable for 10 instruments
- One silicone mat



www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.co.uk



© August 2017 Duckworth & Kent
05.12.17

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.

at the Leading Edge