






Suggested Sub Tenon Set

9-556 Kratz Barraquer Speculum

- 14.5mm open blades
- Angled to rest temporally
- Single piece construction




1-501B DK Westcott Style Stitch Scissors (curved)



- Blunt tips, curved blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 116mm



2-510-1E Moorfields Utility Forceps



- Delicate serrated tips
- Straight shafts
- Flat handle, length 115mm
- Design registration number 004383396-0001



www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.co.uk



© December 2015 Duckworth & Kent
05.12.17

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.

at the Leading Edge