

# PalmScan™ P2000FP

## FastPach PACHYMETER

**FREE  
30-DAY  
trial**



- **Accurate**
- **Fast Captures**
- **Reliable**
- **Low Cost**

High Frequency Probe - 50 MHz

Capable of Measuring diseased cornea

Rechargeable Lithium-Ion Battery

Easy to Operate

Portable



Available in  
Silver  
White  
Orange  
Light Blue

### PACHYMETER SPECIFICATIONS

Accuracy	4 $\mu$ m
Resolution	1 $\mu$ m
Range	250-900 $\mu$ m
Power Dissipation	150 mW
Probe Frequency	50 MHz
Battery Life on a single charge	Over 300 Patients

Toll Free (866) 730-0663  
Direct (818) 222-3310  
Fax (818) 337-1952  
www.micromedinc.com  
info@micromedinc.com

Micro Medical Devices, Inc.  
23945 Calabasas Rd., Ste. 110  
Calabasas, CA 91302



Solutions That Fit